

Blue Diamond[®]
GEL PRODUCTS & BEAN BAG POSITIONERS

DAVID SCOTT
COMPANY

Patient Positioning Solutions & Operating Room Accessories

SINCE 1972



alpinemeditica 

alpinemeditica es: Jesús Del Pino Favila PIFJ731124UW5
Oficinas: Jaime Gomez No. 2802 Int. 2 Col. Paseos Del Sol
Zapopan, Jalisco, México, C.P. 45079
Tels: (+52) 33-15-78-28-02 / 33-15-68-16-02
www.alpinemeditica.com.mx



davidscottco.com



Since 1972, David Scott Company has evolved into one of the medical industry's leading manufacturers & distributors of operating room positioning products. Our patient positioning products are represented in all 50 United States and more than 75 countries internationally.

Our nearly 50 years of industry presence and buying power provides us the ability to supply medical positioning products that best suit our customers' needs and budget. If our products look familiar, don't be surprised as you have probably seen and used them. In fact, our products are private labeled and distributed by some of the medical industries largest O.R. table and positioning manufacturers.

In addition to our distribution efforts, we have developed and manufactured a complete line of positioning products aptly named **Blue Diamond®**. Our **Blue Diamond®** patient positioning product line consists of hundreds of products which include viscoelastic polymer gel positioners, Light Cloud™ reusable-disposable foam head positioners, 3 levels of pressure reducing O.R. table pads, Bean Bag positioners as well as a comprehensive line of O.R. accessories.

Unlike other manufacturers, David Scott Company provides you the ability to order products with no minimum purchase. Purchasing in quantities? Ask for a quantity discount as that will help aid in driving down our already low prices. We have worked with nurses and doctors from all around the world to produce one-of-a-kind custom products when the standard is simply not enough.

This catalog is a representation of some of the more popular products we manufacture and distribute. We have access to thousands of operating room products, so if you do not see what you are looking for please visit our website at www.davidscottco.com or call us toll-free at 800-804-0333 and put our knowledgeable staff to work for you.

DAVID SCOTT
COMPANY


phone 508.875.3333 • toll-free 800.804.0333 • fax 508.875.3375
website www.davidscottco.com • email sales@davidscottco.com



Featured Products



Gel Patient Positioners
See Pages 4-9



Peg Board Lateral Positioning System
See Page 17



Table Pads
See Pages 12 & 13



O.R Side Rail Sockets
See Page 20



Stirrups for Lithotomy
See Pages 21-23



Spine Frames
See Page 25



**Blue Diamond® Gel
Bean Bag Positioners**
See Pages 10 & 11



Urology
See Page 21



Step Stools
See Page 33



Hand Surgery Tables
See Page 16



Roller Boards
See Page 18



**Light Cloud™
Responsible Positioners**
See Page 14

Blue Diamond® dry viscoelastic polymer gel pads and positioners are designed to simulate human adipose tissue composition. This will greatly reduce the risk of decubitus ulcers caused by sheer friction. Independent studies have proven the superiority of gel as compared to foam and other conventional positioning and pressure reducing products. Proudly made in the USA.

THE CONCERN:

Pressure insults and injury.

THE SOLUTION:

Blue Diamond® Gel positioning aids


BLUE DIAMOND® GEL FACTS:

- 100% Latex-Free
- X-Ray Translucent
- Radiolucent
- Antimicrobial
- Repairable
- Antibacterial
- Plasticizer Free
- Silicone Free
- Reusable
- Easily Cleaned
- Will Not Leak, Ooze, or Evaporate
- Contains No Water
- Provides Pressure Reducing Attributes
- Redistributes Weight
- 2-Year Warranty

Gel Head Positioners



BD2150
Adult Head Donut
8"OD x 3"ID x 1.75"H



BD2150-T
Adult Head Donut (Tall)
7.75"OD x 4"ID x 3"H



BD2155
Neonatal Head Donut
3.25"OD x 1.5"ID x .75"H



BD2160
Medium/Pediatric
Head Donut
5.5"OD x 2.25"ID x 1.25"H



BD2165
Neonatal Horseshoe
Head Donut
3.25"OD x 1.25"ID x .75"H



BD2170
Adult Horseshoe
Head Donut
8"OD x 3"ID x 1.75"H




BD2170-T
Adult Horseshoe
Head Donut (Tall)
7.75"OD x 4"ID x 3"H



BD2175
Medium/Pediatric
Horseshoe Head Donut
5.5"OD x 2.25"ID x 1.25"H



BD2180 / BD2185
Head Pillow w/Centering Dish
9"L x 10"W x 1"H (BD2180)
9"L x 10"W x 1.75"H (BD2185)



BD2190
Contoured Head Gel Pad
7"L x 7"W x 3.5"H (outer edge)



BD2191 / BD2192
Richards Headrest
9"L x 8"W x 3"H (BD2191)
9"L x 8"W x 4.5"H (BD2192)



BD2220
Prone Headrest
w/Tube Channels
8"OD x 3"ID x 3.5"H



BD2225
Prone Gel Headrest
w/Rigid Bottom
11"L x 9"W x 6"H

Gel Foot, Heel, Knee & Stirrup Protection

Having the proper padding, cradling, and support is essential for any procedure. Our product line consists of a variety of Blue Diamond® gel pads which include gel pole pads for traditional candy cane stirrups and gel boot pads for lift assisted stirrups. Utilizing foot cradles and heel lifts aids in patient comfort and circulation.



BD2280

Foot, Heel, and Sole Protector w/Velcro®
15"L x 6"W x .25"H



BD2290

Clamshell Heel Pad w/Velcro®
.5" Thick Clamshell
Each




BD2295

Heel Elevators/Lifts
6"L x 2"W x 1.75"H
Pair



BD2300

Small Knee Crutch Pad
9"L x 11"W x 3/8"H



BD2310 / BD2310-5

Medium Knee Crutch Pad
BD2310: 11"L x 17"W x .25"H
BD2310-5: 11"L x 17"W x .5"H



BD2320

Large Knee Crutch Pad
20"L x 20"W x .5"H



BD2325 / BD2326

Gel Stirrup Boot Liners
Bell Shaped
BD2325: Universal Fit
BD2326: Fits Allen® Pal Pro
Pairs




BD2330 / BD2335

Candy Cane Gel Pads
BD2330: 24" Gel Sheath
BD2335: Gel Ankle Protectors
Pairs



BD4502P

Gel Heel Cups
6"L x 3"W x 1.5"H
Pair



BD4502

Gel Heel Cup
6"L x 3"W x 1.5"H
Each

Gel Arm & Elbow Protection

Bony prominences need extra care. If an arm is tucked it becomes vulnerable to ulcerations. We offer a wide variety of product solutions specific to the arm and elbow that will reduce shear, pressure, and the probability of pressure ulcers. We also offer gel for arm shields, armboards, and armboard pads.



BD2230
Small Armboard
Wide Version
15"L x 6"W x .25"H



BD2235 / BD2235-58
Full Size Armboard Pad
w/Rounded End
BD2235: 26"L x 6"W X .25"H
BD2235-58: 26"L x 6"W X 5/8"H



BD2240
Large Full Armboard Pad
24"L x 4.5"W x 3/8"H



BD2241
Lateral Armboard Pad
w/Velcro®
15"L x 6"W x .25"H



BD2245
Toboggan Arm Gel Pad
w/Velcro®
11"L x 17"W x .25"H
(Gel Pad Only)



BD2260 / BD2265
Full Size Elbow Nerve Protector
w/Velcro®
BD2260: 11"L x 17"W x .25"H (Standard)
BD2265: 17"L x 24"W x .5"H (Bariatric)



BD2270
Oval Ulnar Nerve
Protector w/Velcro®
13"L x 6"W x .5"H



BD4307
Small Contoured Armboard
15"L x 5.5"W x 2"H



BD4310
Large Contoured Armboard
20"L x 5.75"W x 2.25"H

Gel Contoured Positioners & Gel Sandbags

Commonly referred to as sandbags, ours are filled with Blue Diamond® Gel which offers pressure reduction compared to the harsh rigidity of an actual sandbag. These positioners are great for tucking and wedging.




BD2375 / BD2370 / BD2360
Gel Filled Sandbag
BD2375: 20"L x 4"W x 2"H (Large)
BD2370: 14"L x 4.5"W x 2.25"H (Med)
BD2360: 10"L x 4.5"W x 2"H (Small)




BD2371
Jumbo Gel Sandbag
13"L x 7"W x 3"H




BD2372
XXL Gel Sandbag
15"L x 8"W x 3"H




BD2380
Sloped Thoracic
Gel Positioner
(From 2"H to 1.5"H)
16"L x 10"W



BD2400
Large Gel Roll Cover
w/Velcro®
20"L x 20"W x .5"H



BD2405
Small Gel Roll Cover
w/Velcro®
17"L x 11"W x .5"H



BD2465
Trapezoid Gel
Positioner
10"L x 7"W x 3.5"H



BD2470
Gel Chest Positioner
Left Side
28"/23L" x 8"W x 2"H



BD2480
Gel Chest Positioner
Right Side
28"/23L" x 8"W x 2"H



BD2485
Pair of Gel Chest Positioners
Left & Right Sides
Same as BD2470 & BD2480

Adult Gel Chest Rolls

Blue Diamond® roll positioners are designed to offer versatility. Our rolls are rounded on the end to allow for axillary usage. The thermoformed skin is very forgiving yet extremely durable. These rolls are perfect for wedging, elevating, distending, padding, and keeping the patient in proper position without fear of flattening the pad or rolling over it.



BD2410
Versatile Roll
12"L x 4"W x 3.75"H



BD2415
Versatile Roll
14"L x 6"W x 3"H
(Dome Shaped)



BD2420
Versatile Roll
16"L x 4"W x 3.75"H



BD2430
Versatile Roll
20"L x 4"W x 3.75"H



BD2440
Versatile Roll
13"L x 7"W x 5.5"H



BD2450
Versatile Roll
16"L x 7"W x 5.5"H



BD2460
Versatile Roll
20"L x 7"W x 5.5"H

Pediatric/Infant Gel Chest Rolls

In addition to our adult gel chest rolls, we offer a complete line of pediatric/infant gel chest roll positioners. These rolls are perfect for wedging, elevating, distending, padding and keeping your pediatric patient in the proper position without fear of flattening the pad or rolling over it.



BD2490
Pediatric Versatile Roll
6"L x 2"W x 2"H



BD2495
Pediatric Versatile Roll
8"L x 2"W x 2"H



BD2510
Pediatric Versatile Roll
10"L x 2"W x 2"H



BD2520
Pediatric Versatile Roll
12"L x 2"W x 2"H



BD2525
Pediatric Versatile Roll
4"L x 3"W x 3"H



BD2530
Pediatric Versatile Roll
6"L x 3"W x 3"H



BD2540
Pediatric Versatile Roll
8"L x 3"W x 3"H




BD2550
Pediatric Versatile Roll
10"L x 3"W x 3"H




BD2560
Pediatric Versatile Roll
12"L x 3"W x 3"H

Gel Operating Table Pads


These Blue Diamond® gel pads are an excellent way to offer extended pressure management especially for cases lasting longer than one hour. Our Blue Diamond® gel table pads cover large surface areas and are used in a wide range of procedures. Can be used in conjunction with foam O.R. table pads on pages 12 & 13.



BD0100-24505
Extra Large & Wide Gel Pad
50"L x 24"W x .5"H
Each



BD2100
3/4 Table Pad
46"L x 20"W x .5"H



BD2120
3/4 Table Pad w/GU Cutout
46"L x 20"W x .5"H



BD2110 / BD2110-25
Full Table Pad
BD2110: 72"L x 20"W x .5"H
BD2110-25: 72"L x 20"W x .25"H



BD2130
3 Piece Table Pad
(2 GU Cutouts)



BD2135
3 Piece Table Pad
(w/GU Cutout on Body Section)




BD2130-1
Head Section of 3 Piece Set
10"L x 20"W x .5"H



BD2130-2
Body Section of 3 Piece Set
(w/GU Cutout)
40"L x 20"W x .5"H



BD2130-3
Foot Section of 3 Piece Set
(w/GU Cutout)
20"L x 20"W x .5"H



BD2140 / BD2140-25
Square Hip/Foot Pad
BD2140: 20"L x 20"W x .5"H
BD2140-25: 20"L x 20"W x .25"H

Gel Specialty Pads

Blue Diamond® pads and positioners can be utilized for many different applications. There are numerous specialty frames, tables, and positioning devices that would benefit from our Blue Diamond® gel pads. If you need help determining the right product please contact us.



BD0100-100705
Gel Pad
10"L x 7"W x .5"H



BD0100-103650
Gel Pad (Backboard Pad)
10"L x 36"W x .5"H




BD0100-172405
Gel Pad
17"L x 24"W x .5"H



BD2250
Shoulder Brace Pads w/Velcro®
7"L x 5"W x .5"H Gel Area
w/ 5" Wraps
Pair



BD2462
Orthopedic Hip Positioner Gel
with Foam Core
16"L x 9"W x 3"H



BD2570
Fracture Table Hip Pad
7"L x 5"W x .5"H



BD2580
Fracture Table Post Pad
w/Velcro® Wraps
7"L x 5"W x .5"H



BD2610
Fracture Table Hip Pad
Accepts up to 2.5" Post
10"L x 8"W x .5"H



BD2640
Wheelchair Seat Cushion
16"L x 18"W x 5/8"H



BD2650 / BD2655
Kambin/Wilson Frame Pads w/Velcro®
BD2650: 29"L x 10"W x .5"H
BD2655: 29"L x 10"W x 1"H
Pair



BD2670
McGuire/IMP Post Caps
3"L x 3"W x .5" Thick
Pair



BD2700
Hip/Sacral Pad Set
1-BD2800 & 1-BD2810
Set



BD2800
McGuire™ / IMP® Base Pad
11"L x 17"W x 1/2"H
Each



BD2810
McGuire™ / IMP® Back Pad
9"L x 11"W x 1/2"H
Each

Custom sizes available, call for prices
800-804-0333

Traditional & Gel Bean Bag Positioners

The safest and most effective way to position patients in a variety of positions such as:

**Lateral • Shoulder Support • Prone • Thyroid • Kidney
Head Immobilization • Frog Leg • Full Length Body Support**


We manufacturer the world’s only pressure reducing Gel Bean Bag positioners on the market (protected by US Patent # 6,226,820).

- Replaceable Valve System
- One-Click Retain/Release Valve
- Valve Safety Clip Included
- Direct Connect to Wall Suction
- Repairable if Punctured or Torn
- Optimal Site Access
- Repels Betadine Staining
- Pressure Insult Prevention
- 100% Latex-Free
- Silicone Free
- Easy to Clean & Disinfect
- Reusable
- Made of Urethane Film (No PVC)
- 6 Sizes Available in Traditional & Gel
- Made in the USA
- 2-Year Warranty



Positioning as simple as 1-2-3-4


1. Place Bean Bag on table, smooth out the beads then place patient in position
2. Mold Bean Bag around patient providing optimal access to surgical site
3. Add suction line until desired firmness to hold the patient in place
4. Simply click our valve, attach the safety clip then remove suction line



BD-BB1020
Traditional Bean Bag
White Valve Replaceable
20"L x 10"W



BD-BB1020-G
Gel Bean Bag
Integral Gel Layer
White Valve Replaceable
20"L x 10"W



BD-BB2020
Traditional Bean Bag
White Valve Replaceable
20"L x 20"W



BD-BB2020-G
Gel Bean Bag
Integral Gel Layer
White Valve Replaceable
20"L x 20"W



BD-BB3040
Traditional Bean Bag
w/Shoulder Cutout
White Valve Replaceable
40"L x 30"W



BD-BB3040-G
Gel Bean Bag
w/Shoulder Cutout
Integral Gel Layer
White Valve Replaceable
40"L x 30"W



BD-BB4036
Traditional Bean Bag
w/Shoulder Cutout
White Valve Replaceable
36"L x 40"W



BD-BB4036-G
Gel Bean Bag
w/Shoulder Cutout
Integral Gel Layer
White Valve Replaceable
36"L x 40"W



BD-BB4646
Traditional Bariatric Bean Bag
w/Shoulder Cutout
White Valve Replaceable
46"L x 46"W



BD-BB4646-G
Gel Bariatric Bean Bag
w/Shoulder Cutout
Integral Gel Layer
White Valve Replaceable
46"L x 46"W




BD-BB3076
Traditional Full Length
Bean Bag
White Valve Replaceable
76"L x 30"W



BD-BB3076-G
Gel Full Length Bean Bag
Integral Gel Layer
White Valve Replaceable
76"L x 30"W

Bean Bag Options & Accessories



BD3500
Gel Overlay w/Cutout for
Bean Bag (BD-BB3040)
40"L x 30"W x .25"H



BD2140-25
Gel Overlay for Bean
Bag (BD-BB2020)
20"L x 20"W x .25"H



BD2110-25
Gel Overlay for Bean
Bag (BD-BB3076)
72"L x 20"W x .25"H



BD-BB176226
White Valve Replacement Kit
w/Safety Clip



BB-CRT
Bean Bag Crimp Tool



BB-HP3
Bean Bag Hand Pump



BD-BB-R
Bean Bag/Gel Bean Bag Repair Patch
3" x 6" (Cut-to-size needed)


SheerCare Series (SG): High Density Foam covered by a durable yet soft Fusion III 4-Way stretch material. Similar to your original pad.

UltraCare Series (UC): Same as SheerCare Series with added top layer of ½" Visco Memory Foam beneath the Fusion III covering. Pressure Management Pad, great for longer duration cases.

Combo Series (C): Same as UltraCare Series with an additional layer of **Blue Diamond®** Gel beneath the covering. Pressure Management Pad, great for longer duration cases.

OR Table Manufacturer/Table #	Thickness	SheerCare	UltraCare	Combo
Berchtold B810, B710	2"	OP-810-SG	OP-810-UC	OP-810-C
	3"	OP-810-SG-3	OP-810-UC-3	OP-810-C-3
	4"	OP-810-SG-4	OP-810-UC-4	OP-810-C-4
Berchtold Operon D850, D750 (long foot)	2"	OP-850L-SG	OP-850L-UC	OP-850L-C
	3"	OP-850L-SG-3	OP-850L-UC-3	OP-850L-C-3
	4"	OP-850L-SG-4	OP-850L-UC-4	OP-850L-C-4
Berchtold Operon D850, D750 (short foot)	2"	OP-850S-SG	OP-850S-UC	OP-850S-C
	3"	OP-850S-SG-3	OP-850S-UC-3	OP-850S-C-3
	4"	OP-850S-SG-4	OP-850S-UC-4	OP-850S-C-4
Maquet 1130	2"	1130-SG	1130-UC	1130-C
	3"	1130-SG-3	1130-UC-3	1130-C-3
	4"	1130-SG-4	1130-UC-4	1130-C-4
Maquet 1130 split leg	2"	1130L-SG	1130L-UC	1130L-C
	3"	1130L-SG-3	1130L-UC-3	1130L-C-3
	4"	1130L-SG-4	1130L-UC-4	1130L-C-4
Maquet 1131.02	2"	1131.02-SG	1131.02-UC	1131.02-C
	3"	1131.02-SG-3	1131.02-UC-3	1131.02-C-3
	4"	1131.02-SG-4	1131.02-UC-4	1131.02-C-4
Maquet 1132 AlphaStar (4 Piece Set)	2"	1132.04-SG	1132.04-UC	1132.04-C
	3"	1132.04-SG-3	1132.04-UC-3	1132.04-C-3
	4"	1132.04-SG-4	1132.04-UC-4	1132.04-C-4
Maquet 1132 AlphaStar (5 Piece Set)	2"	1132-SG	1132-UC	1132-C
	3"	1132-SG-3	1132-UC-3	1132-C-3
	4"	1132-SG-4	1132-UC-4	1132-C-4
Maquet 1133 Alphamaxx (long torso)	2"	1133L-SG	1133L-UC	1133L-C
	3"	1133L-SG-3	1133L-UC-3	1133L-C-3
	4"	1133L-SG-4	1133L-UC-4	1133L-C-4
Maquet 1133 Alphamaxx (short torso)	2"	1133S-SG	1133S-UC	1133S-C
	3"	1133S-SG-3	1133S-UC-3	1133S-C-3
	4"	1133S-SG-4	1133S-UC-4	1133S-C-4
Gettinge/Shampaine 4900B / 4900E	2"	4900-SG	4900-UC	4900-C
	3"	4900-SG-3	4900-UC-3	4900-C-3
	4"	4900-SG-4	4900-UC-4	4900-C-4
Gettinge/Shampaine 5100B / 5100E	2"	5100-SG	5100-UC	5100-C
	3"	5100-SG-3	5100-UC-3	5100-C-3
	4"	5100-SG-4	5100-UC-4	5100-C-4

SheerCare Series (SG)



UltraCare Series (UC)



Combo Series (C)



OR Table Manufacturer/Table #	Thickness	SheerCare	UltraCare	Combo
Skytron 5001/6001/6002 6500/6600/HD	2"	6500-SG	6500-UC	6500-C
	3"	6500-SG-3	6500-UC-3	6500-C-3
	4"	6500-SG-4	6500-UC-4	6500-C-4
Skytron 3100 (2-Piece Set)	2"	3100-SG	3100-UC	3100-C
	3"	3100-SG-3	3100-UC-3	3100-C-3
	4"	3100-SG-4	3100-UC-4	3100-C-4
Skytron 3500/3500B	2"	3500-SG	3500-UC	3500-C
	3"	3500-SG-3	3500-UC-3	3500-C-3
	4"	3500-SG-4	3500-UC-4	3500-C-4
Skytron 3501B/3502 3600B/3602	2"	3501B-SG	3501B-UC	3501B-C
	3"	3501B-SG-3	3501B-UC-3	3501B-C-3
	4"	3501B-SG-4	3501B-UC-4	3501B-C-4
Skytron 3600	2"	3600-SG	3600-UC	3600-C
	3"	3600-SG-3	3600-UC-3	3600-C-3
	4"	3600-SG-4	3600-UC-4	3600-C-4
Skytron 6300/6302	2"	6300-SG	6300-UC	6300-C
	3"	6300-SG-3	6300-UC-3	6300-C-3
	4"	6300-SG-4	6300-UC-4	6300-C-4
Skytron 6700/6701	2"	6700-SG	6700-UC	6700-C
	3"	6700-SG-3	6700-UC-3	6700-C-3
	4"	6700-SG-4	6700-UC-4	6700-C-4
Steris/Amsco 1080, 2080	2"	1300-SG	1300-UC	1300-C
	3"	1300-SG-3	1300-UC-3	1300-C-3
	4"	1300-SG-4	1300-UC-4	1300-C-4
Steris/Amsco 3080,3085	2"	1302-SG	1302-UC	1302-C
	3"	1302-SG-3	1302-UC-3	1302-C-3
	4"	1302-SG-4	1302-UC-4	1302-C-4
Steris C-Max 4085,5085	2"	CMAX-SG	CMAX-UC	CMAX-C
	3"	CMAX-SG-3	CMAX-UC-3	CMAX-C-3
	4"	CMAX-SG-4	CMAX-UC-4	CMAX-C-4
Steris 5085SRT	2"	5085SRT-SG	5085SRT-UC	5085SRT-C
	3"	5085SRT-SG-3	5085SRT-UC-3	5085SRT-C-3
	4"	5085SRT-SG-4	5085SRT-UC-4	5085SRT-C-4
Stryker Vertier	3"	N/A	0788-200-003	N/A

All pads for all tables. Call 800-804-0333 to enquire about CT, MRI, Radiology, Stretcher, and Custom foam pads

Foam Reusable Armboard Pads

All of our foam reusable armboard pads are offered in the same pressure reducing series as our O.R. table pads. Our reusable flexible prone armboard pad is a great alternative to the disposable version. Armboard pads are available in standard, visco memory foam and gel (See top of pages 12 & 13 for a detailed descriptions of our SheerCare, UltraCare and Combo series offerings).



1700-SG
Flexible Contoured
Armboard Positioner
24"L x 6"W
Each



BD-RLP
Reusable Laminectomy Pad
6"L x 6"W x 5"H
Each



BD2010
Gel/Foam Combo Series
Tri-Level Armboard Pad
Contoured
26"L x 6"W
Each



1400-SG
SheerCare Series Armboard Pad
26"L x 6"W
Available in 2", 3", 4", 5", 6" Thickness
Each



1400-UC
UltraCare Series Armboard Pad
26"L x 6"W
Available in 2", 3", 4", 5", 6" Thickness
Each



BD2005
Gel/Foam Combo Series
26"L x 6"W
Available in 2", 3", 4", 5", 6" Thickness
Each

Light Cloud™ Resposable Head/Prone/Heel Positioners

Our Light Cloud™ reusable/disposable foam product line offers an excellent cost saving alternative while providing pressure management attributes. These products are truly unique and wonderful. The doughy soft viscoelastic memory foam is molded to shape and the thin flexible coating is puncture and tear resistant.

**Resposable (Reusable/Disposable) doughy soft, skinned, molded memory foam.
Light Cloud™ positioners can be easily cleaned and reused until the
skin has become compromised.**



LC150
Head Donuts
Overall Dimensions
9"W x 2.5"H
Center Ring 4"
Case of 3



LC250
Prone Positioner
Overall Dimensions
9"L x 9.5"W
Outer Wing Height is 6"H
Chin & Forehead Height is 3.5"H
Each



LC300
Prone Positioners
Case of 3



LC500
Heel Protectors
6"W x 3"L x 1.5"H
Pair

Armboards

All armboards are available with 2", 3", 4", 5", or 6" thick pads. (See top of pages 12 & 13 for descriptions of our SheerCare, UltraCare, or Combo Series foam configurations.)



BDAB26
Ratcheting Phenolic Armboard
21-Positions
26"L x 6"W



BDAB-BARI
Bariatric Armboard
Drop Latch Style
26"L x 10"W



BDAB-CF
Carbon Fiber Underpad
Adjustable Positions
26"L x 6"W



BDAB-LO
Lollipop Underpad Mount
(Pad Sold Separately)
41"L x 26"W x 14"H



BD0102
Gearless Multi-Axis Armboard
Adjustable Height
(Includes 2" Pad)
26"L x 6"W x 2"H



DSC-ABMT
Multi-Task Lateral Armboard
(Includes 2" Pad)
20"L x 6"W x 2"H



DSC-10-380-US
Reison Armboard
Heavy Duty, Stainless Steel
Adjustable, Integral Clamp
(Pad Sold Separately)



DSC-FISSO
Lateral Arm Positioner
3 Rotating Joints
Locks in Place, Smooth Rotation
(Includes Deluxe Pad)



DSC-C-902925
C-Thru Transparent Arm Ladder
4 Adjustable Slots
Radiolucent



DSC-7500-K
Cross Arm Support
Multi-Position Support
Cradle Measures 6"W x 20"L



DSC-7520-K
Cross Cloth Roll
Disposable Sleeves
10 Yard Roll



DSC-800-0296
Armboard Holder
Holds 2 Armboards
Wall Mounted

Arm Shields

Available in Angled and Sled-Style, our arm shields come in standard sizes and can be custom designed to your specifications. Gel protective pads available for most arm shields.



BDAS6
L-Shaped Arm Shield
10"W x 7"D x 6"H
Optional Gel Liner BD2300-AS6
(Pair)



BDAS8
L-Shaped Arm Shield
16"W x 7"D x 8"H
Optional Gel Liner BD2245
(Pair)



BDAS10
L-Shaped Arm Shield
17"W x 7"D x 10"H
Optional Gel Liner BD2245
(Pair)



BD2245
Toboggan Arm Gel Pad
w/Velcro®
11"L x 17"W x .25"H
(Gel Pad Only)



BDAS14125
Sled Style Arm Shield
14"W x 12"D x 5"H
(Pair)



BDAS17
Sled Style Arm Shield
17"W x 10"D x 6"H
(Pair)



BDAS181165
Sled Style Arm Shield
18"W x 11"D x 6.5"H
(Pair)



BDAS221165
Sled Style Arm Shield
22"W x 11"D x 6.5"H
(Pair)



BDAS271165
Sled Style Arm Shield
27"W x 11"D x 6.5"H
(Pair)

Arm & Hand Surgery Tables

We offer both Phenolic and Carbon Fiber hand tables, some have integral clamps and others have a unique under table pad mount which require no clamps at all. Styles include rectangular, hourglass, ironing board, and several tables have the choice of adding an optional side rail. We offer a 2-year warranty against defects in materials and workmanship.



BD301
T-Shaped Phenolic Table
OR Surface 12" x 30"
Includes Clamps & 2" Pad



BD310 (Single Post)
BD311 (T-Post)
Rectangular Phenolic Table
OR Surface 18" x 32"
Includes Clamps & 2" Pad



BD310-MB (Mobile Base)
Free Standing Hand Table
OR Surface 18" x 32"
Height Adjust 30"– 42"
w/Casters & Includes 2" Pad
3 Mounting Options:
Roll-up, Underpad, or Side Rail Mount



BD315 (Single Post)
BD316 (T-Post)
Rectangular Carbon
Fiber Table
OR Surface 18" x 32"
Includes Clamps & 2" Pad



BD320 (Single Post)
BD321 (T-Post)
Hourglass Phenolic Table
OR Surface 15" x 32" x 24" (foot)
Includes Clamps & 2" Pad



BD320-MB (Mobile Base)
Free Standing Hourglass Hand Table
OR Surface 15" x 32"
Height Adjust 30" – 42"
w/Casters & Includes 2" Pad
Roll-up, Underpad, or Side Rail Mount



BD325 (Single Post)
BD326 (T-Post)
Hourglass Carbon Fiber Hand Table
OR Surface 15" x 32" x 24" (foot)
Includes Clamps & 2" Pad



BD330
Rectangle Phenolic
Legless Hand Table
OR Surface 18" x 32"
Includes Clamps & 2" Pad



BD340 (Single Post)
BD341 (T-Post)
Rectangle Phenolic
Under Pad Mount
OR Surface 18" x 32"
Includes Clamps & 2" Pad



DSC-CST-2007
Carter Hand Traction Table
Converts from Stable
Traction to Legless
Height Adjustable 24" to 34"
Includes 1" Thick Pad &
Storage Case



DSC-C-902155
DSC-C-902707
Transparent Under Mount
C-Thru Table
DSC-C-902155: 45"L x 15"W
DSC-C-902707: 45"L x 11"W
Pad Sold Separately

Peg Board & Accessories


Designed for lateral positioning, our peg board system provides surgeons a wide range of positioning options. We offer complete systems as well as individual components.




WSHP0100
Peg Board System
System includes Peg Board,
4-Short & 4 Long Pegs,
2 Safety Clamps &
BD2100 Gel Pad (46"L x 20"W)
(Requires 2 BD-RC Side Rail Clamps)

WSHP0100-M
Peg Board System w/Metal Pegs
System includes Peg Board,
4-Short & 4 Long Metal Pegs,
2 Safety Clamps
BD2100 Gel Pad (46"L x 20"W)
(Requires 2 BD-RC Side Rail Clamps)


WSHP0100-LP
Peg Board System w/Large Gel Pad
System includes Peg Board,
4-Short & 4 Long Pegs,
2 Safety Clamps &
BD0100-24505 Gel Pad (50"L x 24"W)
(Requires 2 BD-RC Side Rail Clamps)



WSHP0101
Peg Board Only
48"L x 19"W x 1"H
Each



BD2100
Gel Pad Overlay
46"L x 20"W x .5"H
Each



BD0100-24505
Extra Large & Wide Gel Pad
50"L x 24"W x .5"H
Each



WSHP0104-PR
Peg Board Safety Clamps
Attaches from peg board to side rail
Pair



WSHP0102
6" Radiolucent Pegs
Set of 4



WSHP0102-M
6" Metal Pegs
Set of 4



WSHP0103
9" Radiolucent Pegs
Set of 4



WSHP0103-M
9" Metal Pegs
Set of 4



WSHP-PE
4" Peg Extender
Each



BD0100-2580
Gel Wrap for 6" Pegs
Each



BD2300-8P
Gel Wrap for 9" Pegs
Each



BD2300-10P
Gel Wrap for 10" Pegs
Each



BD2310-5-12P
Gel Wrap for 12" Pegs
Each



BD-RC
Peg Board Clamp
Blade Style
Each

Patient Roller Boards

There is no substitute. While many try to knockoff, we distribute the original Samarit roller board. No shortcuts in materials and workmanship, only the best for over 30 years.



Samarit SMART Roll Board Facts:

- Flame Retardant Surface
- Anti-Bacterial
- Latex-Free / Skin Friendly
- Requires No Lifting
- Easy Patient Transfer
- 3-Year Warranty



DSC-RB4100
Standard Smart Roller Board
70"L x 20"W Foldable
6 lbs (2.7 kg)



DSC-RB4500
MRI/CT Roller Board
70"L x 16"W Foldable
5.5 lbs (2.5 kg)



DSC-RB4400
Surgi-Long Roller Board
48"L x 16"W Non-Foldable
3.25 lbs (1.5 kg)



DSC-RB4460
Surgi-Wide Roller Board
48"L x 20"W Non-Foldable
3.25 lbs (1.5 kg)



DSC-RB4450
Surgi-Short Roller Board
40"L x 16"W Non-Foldable
2.75 lbs (1.3 kg)



DSC-RB4300
Mini Long Roller Board
40"L x 20"W Non-Foldable
3.25 lbs (1.5 kg)



DSC-RB4350
Mini Short Roller Board
36"L x 20"W Non-Foldable
3 lbs. (1.4 kg)




DSC-RBR
Roller Board Mounting
Brackets Powder Coated
Includes 1 Top & 1 Bottom



DSC-RB2515 (Small)
DSC-RB3014 (Medium)
DSC-RB6714 (Large)
Roller Boards w/Cover
DSC-RB2515: 25"L x 15"W
DSC-RB3014: 30"L x 14"W
DSC-RB6714: 67"L x 14"W

Transfer Boards



4777XL
Full Length Patient Shifter
72"L x 22"W x 3/16"Thick



4755XL
Standard Patient Shifter
33"L x 22"W x 5/8"Thick



4777-HANDLE
Handles for 4755 & 4777 Shifters
Pair



DSC-WHPS
Universal Wall Hanger
for Patient Shifters
(4755XL & 4777XL)
Adjustable

Patient Restraints

Our patient restraints are available in different lengths and styles. Our catalog represents the restraints that we have found to be most popular. We recommend the addition of a removable **Blue Diamond®** gel layer, which will reduce the risk of shear and pressure points. Custom lengths and styles are available.



SWQTR5X30GY

Disposable Velcro® Hook & Loop
Self-Wrapping
5" x 30 Feet (1 Roll)
See Website for Additional Sizes



SWQTR5X30GY-10

Disposable Velcro® Hook & Loop
Self-Wrapping
5" x 30 Feet (Case of 10 Rolls)
See Website for Additional Sizes



BD700/BD700L/BD700XL

Simple Rubber Strap
1 Airplane Buckle
BD700: 2" X 72"
BD700L: 2" x 120"
BD700XL: 2" x 132"
Custom Lengths Available



BD800/BD800L/BD800XL

Conductive Patient Restraint
2 Airplane Buckles/All Rubber Mid Panel
BD800: 96" w/18" x 5" Mid Panel
BD800L: 120" w/18" x 5" Mid Panel
BD800XL: 132" w/18" x 5" Mid Panel
Custom Lengths Available




BD2620/BD2620L/BD2620XL

Same as BD800 with Gel Pad
Gel Pad Attaches to Mid Panel
BD2620: 96"
BD2620L: 120"
BD2620XL: 132"
Custom Lengths Available



BD800-PAD/BD800L-PAD/BD800XL-PAD

Same as BD800 w/ Padded Memory Foam Mid Panel
BD800-PAD: 96"
BD800L-PAD: 120"
BD800XL-PAD: 132"



BD8000

Bariatric Double Restraint
Padded 30" x 12" Mid Panel
and w/2 80" Straps



BD810

Hooks for Restraint Straps
Sold as Each



BD800-NM

Non-Metal Version of BD800
96" w/ Velcro® Closures



BD805

Same as BD800
96" w/ Nylon Straps




BD805-NMP

Same as BD805
96" w/ No Mid Panel
2 Airplane Buckles



BD805L-NMP

Same as BD805
120" w/ No Mid Panel
2 Airplane Buckles



DSC-DV12

DaVinci® Gel Restraining Bump
21"L x 2.5"H w/ Gel Bump
(Includes Integral Gel Pad)



DSC-NTS

Thigh Retractors
(Nissen Style)
Pair



DSC-NTS-B

Bariatric Thigh Retractors
(Nissen Style)
Pair

O.R. Side Rails Sockets/Clamps



BD-CS
Stainless Steel Clark Socket
Accepts Up To 5/8"
Round Post



DSC-VCS
Universal Clark Socket
Accepts Both Round & Blade Posts
Attaches Anywhere Along Rail



BD-URC
Universal Rail Clamp
Accepts Both Round & Blade Posts
Attaches Anywhere Along Rail



BD-BC
Stainless Steel
Accepts Blade Post Only



BD-RC
Blade Style Rail Clamp
Accepts Blade Post Only



BD-USRC
Stainless Steel
Drops Anywhere on Rail
Accepts Round Post



DSC-800-0134
Schure Socket XL
Accepts Both Round & Blade Posts
Attaches Anywhere Along Rail



DSC-800-0271
CamLoc Clamp
Accepts Both Round & Blade Posts
Attaches Anywhere Along Rail



BD168-6
6" Add-A-Rail
For Tables w/o Side Rails




DSC-SRE
Side Rail Extender
Adds 8" Side Rail




DSC-800-0338
SpringLoc Clamp
Attaches Anywhere Along Rail
Accepts Blade Posts


Foot Board & Table Width Extenders



BD080
Foot Board Extender
Adjustable Width 18.75"-23.75"
Requires 2 BD-RC Clamps



DSC-TWE
Table Width Extender
4" Extender w/locking integral clamps
DSC-TWE-2 (Includes 2" pad)
DSC-TWE-3 (Includes 3" pad)
DSC-TWE-4 (Includes 4" pad)



DSC-TWE-XL
Bariatric Table Width Extender
8" Extender w/locking integral clamps
DSC-TWEXL-2 (Includes 2" pad)
DSC-TWEXL-3 (Includes 3" pad)
DSC-TWEXL-4 (Includes 4" pad)

Urology Frame & Drapes



DSC-UC
Flexible Urocatch Frame
w/Integral Clamps



DSC-UC31
Sterile Urocatch Drapes
Case of 10



DSC-UC32A
Non-Sterile Urocatch Drapes
Case of 10



DSC-15LC-P
15 ltr Collection Container
with Powder
Case of 10

Lift-Assisted Stirrups



BD100

Lift-Assisted Stirrups
Travels 120°, from +90° to -30°
Abduction & Adduction Done at Handle
Includes Boot Liners & Side Rail Clamps
500 lb. Weight Capacity
Pair



DSC-800-0342-PL

Great White Platinum Stirrups
Range of Motion: -30° to +80°
Abduction +25° & Adduction -9°
Includes Premium Boot Liners
600 lb. Weight Capacity
Optional: Clamps Sold Separately
Pair



DSC-800-0342-PR (No Fin)

Great White Premium Stirrups
Range of Motion: -30° to +80°
Abduction +25° & Adduction -9°
Includes Premium Boot Liners
400 lb. Weight Capacity
Optional: Clamps Sold Separately
Pair



DSC-800-0342-R

Great White Robotic Stirrups
Range of Motion: 65° low to 45° High
Abduction +25° & Adduction -9°
Includes Premium Boot Liners
600 lb. Weight Capacity
Optional: Clamps Sold Separately
Pair



KTEK3100 (With Fin) / KTEK5100 (With Fin)

Lift-Assisted Stirrups
No Velcro® brand fasteners
Range of Motion: -55° to +85°
Abduction +25° & Adduction -9°
KTEK3100: 350 lb. / KTEK5100: 500 lb.
Includes Boot Liners
Optional: Clamps Sold Separately
Pair

Bariatric Lift-Assisted Stirrups



BD150

Bariatric Lithotomy Stirrups
Range of Motion: -33° to +84°
Abduction +25° & Adduction -9°
Includes Premium Boot Liners
800 lb. Weight Capacity
Optional: Clamps Sold Separately
Pair



DSC-800-0364 (With Fin)

World's First Interchangeable Boot System
Includes 3 Sets of Boots w/Pads (Premium, Platinum, & Bariatric)
Range of Motion: -40° to +75°
Abduction +25° & Adduction -9°
1000 lb. Weight Capacity
Optional: Clamps Sold Separately
Pair

Pediatric & Infant Stirrups



DSC-800-0342-PUP (160 lbs.)
DSC-800-0342-KIDS (250 lbs.)

Pediatric E-Z Lift Stirrups
Range of Motion: -35° to +84°
Abduction +25° & Adduction -9°
Includes Boot Liners
Pair



BD010

Infant/Pediatric Stirrups
Flexible for Infinite Positioning
Includes Integral Clamps
Disposable Style Loop Straps (Sold Separately)
Pair



BD010-S

Infant Stirrups Straps
For BD010 Stirrups
Case of 10

Basic Stirrups



BD007

Basic Stirrup Leg Holder
Upgrade from Candy Canes
Requires 2 Clark Sockets (BD-CS)
Pair



BD-CC

Candy Cane Stirrups
Adjusts up to 40" off the Table
Includes Cotton Foot Straps
Requires 2 Clark Sockets (BD-CS)
Optional: Gel Pads Sheaths (BD2330)
Gel Ankle Straps (BD2335)
Pair



BD2330 / BD2335

Candy Cane Gel Pads
BD2330: Gel Sheath
BD2335: Gel Ankle Protectors
Both Items Sold as Pairs

Stirrup Accessories



BB-UFBP-4962
Replacement Liners
UltraFin® Stirrups
Pair



BD-YFBP
Replacement Liners
(Reinforced at Seams)
Yellofin™ Stirrups
Pair



BD100-L
Replacement Liners
BD100 Stirrups
Pair



BD150-L
Replacement Liners
BD150 Stirrups
Pair



BD-AFHBP
Replacement Liners
Allen® Style Foot Huggers
Pair



BD2325
Gel Stirrup Liner
Universal
Pair



BD2326
Gel Stirrup Liner
Allen® Pal Stirrups
Pair



BD2330
Gel Candy Cane Sheaths
Candy Cane Stirrups
Pair



BD2335
Gel Ankle Straps
Candy Cane Ankle Straps
Pair



BD-CC-H
Candy Cane Hook
(Hook Only)
Each



BD-CC-ST
Candy Cane Nylon Straps
(Nylon Strap Only)
Each



DSC-SC
Stirrup Cart
Stainless Steel
Holds 2 Pairs



DSC-800-0075
Stirrup Cart Holder
DSC-800-0342's & DSC-800-0364
See Pages 21 & 22



BD010-S
Infant Stirrups Straps
For BD010 Stirrups
Case of 10



DSC-800-0284
Stirrup Wall Rack
Each

Stirrup Repairs

We offer repair services for all stirrups at a fraction of the cost you would pay from the original manufacturers. For a fixed price, we will thoroughly examine your stirrup fixing what is necessary to bring it back to full functionality.

**Lift-Assisted Gas Cylinder • Handle Repairs • Boot Repairs
Lockdown Repair/Replace • Pivot Joints • Clean/Adjust • Velcro®**



Shoulder Positioners



BD-SS
Shoulder Braces
Adjustable on Linear Plane
Integral Clamps



DSC-800-0100
Shoulder Supports
Adjusts Vertically
(Requires 2 Clark Sockets)



DSC-800-0280
XPS Surgical Shoulder Positioner
Similar in Functionality to Arthrex Trimano®
Smooth Effortless Positioning
(See Website for More Details)



DSC-800-0004
Powered Shoulder Beach Chair
Only Self-Powered Beach Chair on Market
Articulation Range 0-90° & adjusts in 8 seconds
Includes 5-piece Deluxe 2" Pad Set
500 lb. Weight Capacity
(See Website for More Details)



DSC-800-0142
E-Z Lift Shoulder Beach Chair
Powered Piston System for Easy Adjust
Articulation Range 0-90°
Includes 5-piece Deluxe 2" Pad Set
500 lb. Weight Capacity
(See Website for More Details)



BD050-HCS
Head & Chin Straps
(Compatible w/Skytron®, OSI®,
Steris®, Berchtold® & Schuremed™)
Case of 10

Lateral Positioners




BD-BB3040-G
Gel Bean Bag
w/shoulder cutout
Integral Gel layer
40"L x 30"W
(See Pages 10 & 11)



DSC-800-0169
Lateral Shoulder Holder
Easy Control/Lockdown for
Internal & External Rotation
Apply Tension w/Knob from 0-15 lbs.
(No Sandbags Needed)



DSC-UPB-DLX
3-Pc lateral Brace System
System Includes Clamps
Adjusts Vertically, Horizontally &
Diagonally



WSHP0100
Peg Board System
System includes Peg Board,
4-Short & 4-Long Pegs
2 Safety Clamps & BD2100 Gel Pad

Hip Positioners



DSC-4150-00
Stulberg Hip Positioner
Attaches Directly to OR Table
(See Website for More Details)



DSC-4050
Wixson Hip Positioner
Total Hip Positioner
(See Website for More Details)


Spinal Frames



WF-1000

Wilson Frame (Radiolucent)

Refurbished to Brand New Condition
Includes Foam Pads & Transport Case
300 lb. Weight Capacity
(See Website for More Details)



DSC-800-0235

Schure Spine Frame

Easily Attached to OR Table
Includes Foam Pads
500 lb. Weight Capacity
(See Website for More Details)



DSC-800-0074-SF

Spine Frame Dolly


Stainless Steel, Easy Transport



DSC-CS-2608

Bow® & Wilson Frame Covers

Waterproof Foam Sleeves
Disposable, Case of 5



BD2650/BD2655

Kambin/Wilson Frame Gel Pads w/Velcro®

BD2650: 29"L x 10"W x .5"H
BD2655: 29"L x 10"W x 1"H

Repair Patches



BDR001

Gel Repair Patch

6"L x 3"W
Cut To Size Needed
(For all gel pads pages 4-9)



BDR002

Gel Repair Patch

6"L x 24"W
Cut To Size Needed
(For all gel pads pages 4-9)



BD-BB-R

Bean Bag/Gel Bean Bag Repair Patch

3"L x 6"W
Cut To Size Needed
(For all bean bags pages 10 & 11)

Velcro®



VEL-KIT

Adhesive Backed Velcro®
12"L x 2"W Strips
Includes 2-Hook & 2-Loop



191195 / 191245

Adhesive Backed Velcro®
25 Yard Roll – 2" Width
191195: 25 Yards Loop
191245: 25 Yards Hook



SWQTR5X30GY

Disposable Velcro® Hook & Loop
Self-Wrapping
5" x 30 Feet (1Roll)
See Website for Additional Sizes



SWQTR5X30GY-10

Disposable Velcro® Hook & Loop
Self-Wrapping
5" x 30 Feet (Case of 10 Rolls)
See Website for Additional Sizes

Anesthesia Screens




BD200

Malleable Anesthesia Screen
Integral Side Rail Clamps



BD200-NC

Malleable Anesthesia Screen
5/8" Round Posts
No Clamps
(Requires 2 BD-CS Clamps)



BD202

Flexible Table Mount
Integral Clamps
Ideal For Imaging Tables w/No Side Rails



BD205

Standard Rigid Anesthesia Screen
5/8" Round (Requires 1 BD-CS Clamp)



BD206

Adjustable Height Anesthesia Screen
5/8" Round (Requires 1 BD-CS Clamp)

Nerve Stimulators



DSC-SSP

SunStim™ Pro
Nerve Stimulator
See Website for Details



DSC-SSP-BANDS

Color bands for the DSC-SSP
See Website for Color Options



10023

Microstim Nerve Stimulator
See Website for Details



7100

Microstim Nerve Plus Stimulator
See Website for Details



100057

6' Alligator Leads
For Nerve Stimulators



100067

6' Alligator Leads w/Snap
Connectors
For Nerve Stimulators



100262

Diagnostic Ball Probes
For Nerve Stimulators

Antifog



5010

Blue Diamond® Antifog Kits
Adhesive Sponge and 6gm. Vial
20 Sterile Units

Stick-On Thermometer



5030

Stick-On Thermometers
Package of 100

Knee & Leg Positioners with Optional Accessories



BD170
Knee Crutch Supports
(Requires 2 Clark Sockets BD-CS)
Pair



DSC-2760-00
Fromm Knee Triangles
(XS Shown is Sold Separately)
Set of 3



DSC-TKS-R
Total Knee Stabilizer
(Requires Clark Sockets BD-CS)
Each



DSC-TKS-PAD
Foam Replacement Pads
Case of 10




DSC-800-0141 (24" base)
DSC-800-0171 (30" base)
Total Knee Positioner
Smooth Motion, Locks in Position
Autoclavable



DSC-2630
Robb Leg/Knee Positioner
Steam & Gas Sterilizable
Aluminum Footplate
Includes 3 Disposable Sterile Pads



DSC-112001
Hermanator Foam Knee Pad
Sterile/Case of 10
(Works with Alvarado, Robb Leg,
DeMayo™ and DSC-800-0141/0171)




BD2310/BD2310-5
Gel Knee Crutch Pads
11" x 17"
Each



DSC-2620
Stulberg Leg/Knee Positioner
Aluminum Footplate
Includes 3 Disposable
Sterile Pads



DSC-SP
Stress Post & Pad
(Requires 1 Clark Socket)



BD250
Arthroscopic Leg Holder
22 Adjustable Width Positions
Requires 1 Clark Socket
(BD-CS)



Options for BD250:
BD250-P12 Foam Inserts (Case of 12)
BD250-SSP20 Rigid Foam Inserts (Case of 20)
BD2500 Gel Insert for Leg Holder (Each)



BB-WLH
Well Leg Holder
Foam/Reusable
Each



BB-EPP
Extremity Positioner
Foam/Reusable
Each



DSC-FROGGER
Vein Harvesting Positioner
Foam/Reusable
Each



DSC-800-0161
Universal Leg Prepper
(Requires 1 Clark Socket BD-CS)
Each



DSC-SLPS
Single Loop Leg Prepper
Includes Pad
(Requires 1 Clark Socket
BD-CS)



DSC-SLPS-FOAM-PAD
Single Loop Prepper Pads
Wipe-Clean Molded Pads
Case of 12



DSC-610
Double Picket Fence Leg Holder
Preps both legs at once
Adjustable Height up to 30"
Stainless Steel w/Rubber Coated Top
Requires 2 Blade Clamps (BD-BC)



DSC-611
Single Picket Fence Leg Holder
Preps both legs at once
Adjustable Height up to 30"
Stainless Steel w/Rubber Coated Top
Requires 1 Clark Socket (BD-CS)

Bariatric Products

David Scott Company offers a wide variety of bariatric positioning aids. Included in this line are table width extenders, restraint straps, bean bag positioners, lift-assisted stirrups, replacement O.R. table pads, stretcher pads, and our line of **Blue Diamond®** gel positioning products.



BD150

Bariatric Lift Assisted Stirrup w/ Pads
Bariatric Lithotomy Stirrups
Range of Motion: -33° to +84
Abduction +25° & Adduction -9°
Includes Premium Boot Liners
800 lb. Weight Capacity
Optional: Clamps Sold Separately
(See Page 22)




BD8000

Bariatric Double Restraint
w/2 80" Reinforced Straps
Padded 30"L x 12"W Mid Panel
(See Page 19)



BDAB-BARI

Bariatric Armboard
Drop Latch Style
26"L x 10"W
(See Page 15)



DSC-TWEXL-4

8" Table Width Extender w/ 2 Locking Integral Clamps
8"L x 20"W w/ Side Rail & Pad
(See Page 20)



DSC-EJMS

Jumbo Stainless Steel Mayo
Foot Operated Hydraulic Foot Pedal
Height Adjustable 36"-53"
Tray Size 28"L x 20"W
Weight Capacity up to 150 lbs.
(See Page 30)



BD-BB4646-G

Gel Bariatric Bean Bag
w/Shoulder Cutout
(Valve is Replaceable)
46"L x 46"W
(See Pages 10 & 11)



DSC-NTS-B

Bariatric Thigh Retractors
(Nissen Style)
Pair
(See Page 19)



DSC-800-0364 (With Fin)

World's First Interchangeable Boot System
Includes 3 Sets of Boots w/Pads (Premium, Platinum, & Bariatric)
Range of Motion: -40° to +75
Abduction +25° & Adduction -9°
1000 lb. Weight Capacity
Optional: Clamps Sold Separately
(See Page 22)



Blue Diamond® Gel Products

Over 150 Styles
(See Pages 4-9)

Disposable Products

The following is a listing of our most popular disposable products.



DSC-UC31/DSC-UC32a
Urology Disposable Drapes
DSC-UC31: Sterile Drapes
DSC-UC32a: Non-Sterile Drapes
Case of 10



DSC-CS-2600
Comfort Covers
For Allen® Advance Table
6 Sets per Box



DSC-CS-2635/DSC-CS-2639
Comfort Plus/Comfort Ultra
Cushion Set
Allen® Compatible
5 Sets per Box



DSC-CS-2610
Jackson Table Covers
Padded Chest Piece
6 Sets per Case



DSC-CS-2520
Prone One Disposable
Face Cushion w/Mirror
Can be used on all OR Tables
8 Units per Box



DSC-CS-2512 (Adult)
DSC-CS-2514 (Child)
Prone Plus Foam Face Cushion
Compatible w/Proneview®
Case of 10



DSC-5348
OSI® Schlein® Shoulder
Chair Disposable Pad Set
Case of 3



DSC-CS-2535 (7")
DSC-CS-2550 (5")
Prone Face Block Pillows
5" (Case of 12)
7" (Case of 8)



SWQTR
Disposable Velcro® Hook & Loop
Self-Wrapping
See Website for Additional Sizes



DSC-760/DSC-761
Sterile Alumi Hands Immobilizer
DSC-760 (L) 9.5" x 8.5" (8 per Case)
DSC-761 (XL) 13"x10.5" (8 per Case)



LC250/LC300
Light Cloud™ Resposable
Prone Reusable/Disposable
LC250: Each
LC300: Case of 3



LC150
Light Cloud™ Resposable
Head Donut
Reusable/Disposable
Case of 3



BDF-2000/BDF-2000-P
Arm Elevation Pillow
BDF-2000 (Adult) 16"x13.5"x12"
BDF-2000-P (Pedi) 9"x7.75"x4.5"
Pair



DSC-112001
Hermanator Knee Pad
Works w/ Alvarado, Robb Leg
803/803S Universal Fit
Sterile / Case of 10



BD250-P12
Arthroscopic Leg Holder
Foam Pad
Adjustable Closure
Case of 12



LC500
Light Cloud™ Resposable
Heel Protectors
Reusable/Disposable
Pair



DSC-STSNMED1
Safe-T-Secure Trendelenburg
Case of 5
See Website for Details



DSC-VISI
VisiWrap Face Protection
360° Protection (Supine/Beach Chair)
Case of 10



5030
Anesthesia Temp Dots
92°-106°F / 34°-41°C
Liquid Crystal Display
Package of 100



BD050-HCS
Head & Chin Straps
(Compatible w/Skytron®, OSI®,
Steris®, Berchtold® & Schuremed™)
Case of 10

See website for additional disposable products

Mayo Stands



8740SS

OR Table Mount Stainless Steel Mayo
Tray Size 13" x 19"
Height Adjustable to 26"



8750SS/8751SS/8752SS

Double Post Stainless Steel Mayo
Foot Pedal Raises & Lowers Tray
Slow Descend for Safety
Height Adjustable 36"-62"
8750SS: Tray Size 13" x 19"
8751SS: Tray Size 16" x 21"
8752SS: Tray Size 20" x 25"



8867SS/8869SS/8871SS

Double Post Stainless Steel Mayo
Foot Pedal Raises & Lowers Tray
Height Adjustable 39.5"-64"
8867SS: Tray Size 13" x 19"
8869SS: Tray Size 16" x 21"
8871SS: Tray Size 20" x 25"



DSC-730SS/DSC-731SS

Knob Control Stainless Steel Mayo
Low Profile, Easy Height Adjustment
Height Adjustable 36"-62"
DSC-730SS: Tray Size 12-5/8"L x 19"W
DSC-731SS: Tray Size 16.25"L x 21.25"W



DSC-EJMS

Jumbo Stainless Steel Mayo
Foot Operated Hydraulic Foot Pedal
Height Adjustable 36"-53"
Tray Size 28"L x 20"W
Weight Capacity up to 150 lbs.



DSC-OSMS

Oversized Mayo Stand
Turn Crank
Height Adjustable 36"-56"
Tray Size 30"L x 26"W
Weight Capacity up to 100 lbs.

Solution Stands



7087SS-HB

Single Solution Stand
8.5 Quart basin
(No Shelf)



7087SS

Single Solution Stand
8.5 Quart basin
(With Lower Shelf)



8087SS-HB

Double Solution Stand
2 - 8.5 Quart basin
(No Shelf)



8087SS

Double Solution Stand
2 - 8.5 Quart basin
(With Lower Shelf)

All of our durable stainless steel tables (Utility, Instrument, and Nested) are constructed of 16 gauge, 304 Stainless Steel with smooth edges. 3" swivel casters are standard and locking casters are available upon request.

Utility Tables



BD520

Stainless Utility Table
16"W x 20"D x 34"H
1 Drawer w/Shelf



BD521

Stainless Utility Table
16"W x 20"D x 34"H
2 Drawer w/Shelf



BD522

Stainless Utility Table
16"W x 20"D x 34"H
3 Drawer w/Shelf



BD523

Stainless Utility Table
16"W x 20"D x 34"H
4 Drawer w/Shelf

Stainless Steel Instrument Tables



FULL SHELF

SST1620-S	16"W x 20"L x 34"H
SST1630-S	16"W x 30"L x 34"H
SST1833-S	18"W x 33"L x 34"H
SST2036-S	20"W x 36"L x 34"H
SST2048-S	20"W x 48"L x 34"H
SST2436-S	24"W x 36"L x 34"H
SST2448-S	24"W x 48"L x 34"H
SST2460-S	24"W x 60"L x 34"H
SST2472-S	24"W x 72"L x 34"H



H-BRACE

SST1620	16"W x 20"L x 34"H
SST1630	16"W x 30"L x 34"H
SST1833	18"W x 33"L x 34"H
SST2036	20"W x 36"L x 34"H
SST2048	20"W x 48"L x 34"H
SST2436	24"W x 36"L x 34"H
SST2448	24"W x 48"L x 34"H
SST2460	24"W x 60"L x 34"H
SST2472	24"W x 72"L x 34"H

Nested Instrument Tables



SST-N-6	Set of 6 Nested Tables
SST142832-N	14"W x 28"L x 32"H
SST163234-N	16"W x 32"L x 34"H
SST183636-N	18"W x 36"L x 36"H
SST204038-N	20"W x 40"L x 38"H
SST224440-N	22"W x 44"L x 40"H
SST244842-N	24"W x 48"L x 42"H



3" Swivel Casters Come Standard on All Tables
Locking Casters Available Upon Request

Kick Bucket & Hampers

Other styles are available



KB900
Kick Bucket & Frame
12 Qt.



KB920
Kick Bucket & Frame
8 Qt.



CB1710
Square Hamper Stand w/Lid
(Supports 36-42 Gal Bags)
Foot Pedal Operated



BL2018/BL2025
Round Hamper
BL2018: 18" x 35"
BL2025: 25" x 35"



BL2230
Folding Adjustable
Hamper
(Fits 18", 22" & 25" Bags)

IV Poles

IV Poles are available in a variety of configurations (stainless steel, chrome, weighted base, spider base, 2-4-6-8 hook and 2" or 3" casters). Below shows the most popular, please visit davidscottco.com for many more options, configurations, and styles.

Other styles are available



1415SS-4
Stainless Steel IV
5-Leg, 4-Hook
Adjustable 52"-94"



DSC015-4
Stainless Steel IV
5-Leg Powder Coated, 4-Hook
Adjustable 52"-94"
Colors: **Black**, **Charcoal**, **Red**,
Blue, **Yellow**, & **Green**



DSC202-4
Chrome IV Pole
5-Leg, 4-Hook
Adjustable 52"-94"



DSC020-SS-4
DSC020-SS-6
Stainless Steel IV
6-Leg Spider Base
Adjustable 50"-92"
Colors: **Black**, **Charcoal**, **Red**,
Blue, **Yellow**, & **Green**



DSC270-SS-4
Stainless Steel IV
5-Leg Spider Base, 4-Hook
Thumb knob, Adjustable 50"-92"
Colors: **Black**, **Charcoal**, **Red**,
Blue, **Yellow**, & **Green**



BLITZ-4-A
Irrigation IV Tower
4-Channels
Holds up to 8.5 ltr IV Bag
12.75"W x 12.75"D x 77.75"H




1450-HC
Irrigation Tower, Heavy Duty
Available in 4, 8, or 28 Hook
Easy Lift Hand Crank
6-Leg Spider Base



DSC232
IV Pole Power Strip
6 Outlets
Integral Clamp, Direct Mount
Universal Mount

Foot Stools



2678SS
EZ Stacking Step Stool
Stainless Steel
12"L x 18"W x 5"H
Each



2678SS-D
Universal Dolly for Stacking Stools
13"L x 23"W x 39"H
Each



DSC-1170
Orange Step Stools
13"L x 17"W x 5.25"H
Includes 2 Tops & 2 Bottoms/Case



DSC-36361
Ergonomic & Stackable Step Stool
15"L x 11"W x 3"H
Case of 2



DSC-36362
Stainless & Stackable Step Stool
17"L x 11"W x 5"H
Case of 2



DSC-36363
Gray Aluminum Stackable Step Stool (300lb)
17"L x 13"W x 4.625"H
Anti-microbial Coating
Case of 4


Anti-Fatigue Comfort Mats

The ultimate anti-fatigue solution in operating rooms. Our comfort mats have non-slip properties, are available in a variety of colors and sizes and can be custom made to meet your specific needs. Beveled edges, antimicrobial, easy to clean, and impervious to all chemicals are just a few of the benefits our comfort mats provide.




ORTHOMAT

OM-203334	20"W x 33"L x .75"H
OM-202434	20"W x 24"L x .75"H
OM-186034	18"W x 60"L x .75"H
OM-224834	22"W x 48"L x .75"H
OM-227234	22"W x 72"L x .75"H
OM-243634	24"W x 36"L x .75"H
OM-244234	24"W x 42"L x .75"H
OM-363634	36"W x 36"L x .75"H
OM-364234	36"W x 42"L x .75"H
OM-366034	36"W x 60"L x .75"H
OM-404834	40"W x 48"L x .75"H
OM-427234	42"W x 72"L x .75"H
OM-2412034	24"W x 120"L x .75"H



ERGOMAT

EM-3020HE	24"W x 36"L x 7/16"H	Autoclavable
EM-3040HE	36"W x 48"L x 7/16"H	
EM-3060HE	36"W x 60"L x 7/16"H	



Ergomats are available in **Red, Blue, Gray, & Green**

**ANTI-FATIGUE COMFORT MATS ARE AVAILABLE IN CUSTOM SIZES
VISIT OUR WEBSITE FOR MORE SIZES, COLOR OPTIONS & STYLES**

Adjustable Stools



ST411
Chrome, 4-Leg
Adjustable Height 17"-27"
15" Seat, w/Foot Ring



ST1113
ST1113-WB (w/Back)
Chrome, 5-Leg Exam Stool
Swivel Adjustable Height 20"-27"
16" Seat, With or Without Back



ST205
ST206 (w/Back)
Chrome, 5-Leg, Pneumatic Foot
Adjustable Height 18"-24"
16" Seat, With or Without Back



ST1210
ST1212-WB (w/Back)
Chrome, 5-Leg, Pneumatic
Adjustable Height 22"-29"
w/Simple touch of the lever
16" Seat, With or Without Back
(ST1212-WB has 17" Foot Ring)



ST21340B
5-Leg, Pneumatic
Adjustable Height 19"-24"
18.5" Tear Drop Seat & Backrest



ST21435 (w/Arms)
SST21435B (w/o Arms)
5-Leg, Hand Operated Task Stool
Adjustable Height 18.5"-24.75"
18.5" Tear Drop Seat & Backrest



7714-PSS
Stainless Steel Stool
Adjustable Height 19"-31"
14" Padded Seat with Heel Rest



7745SS
Stainless Steel Stool
Adjustable Swivel Height 19"-31"
15" Seat



DSC-15511
DSC-15512 (w/Back)
5-Leg, Hydraulic Stool
Adjustable Height 22.5"-29.5"
15" Seat

Reusable Nylon Cart Covers

All cart covers come with Velcro® closure and squared corners. The standard color is Royal Blue unless specified otherwise.



Royal Blue



Reusable Cart Covers

4001-184854-400D	18"W x 48"L x 54"H
4001-184862-400D	18"W x 48"L x 62"H
4001-184872-400D	18"W x 48"L x 72"H
4001-244854-400D	24"W x 48"L x 54"H
4001-244862-400D	24"W x 48"L x 62"H
4001-244872-400D	24"W x 48"L x 72"H
4001-246054-400D	24"W x 60"L x 54"H
4001-246062-400D	24"W x 60"L x 62"H
4001-246072-400D	24"W x 60"L x 72"H

Cart Covers


See website for additional options such as clear routing pocket, handle holes and clear front flap.

Available in multiple colors:

Royal Blue, Navy, Red, White, Yellow, Kelly & Frost Green

Plastic Sterilization Trays

- Made of Ultem (A high impact resin)
- Trays are stackable
- Pin mats and posts available
- Sterilized by EO, Autoclave, Dry Heat, Cold Solutions
- See through lids
- Multiple sizes available



Plastic Sterilization Trays

6100	Micro tray w/ base, lid, & mat	6" x 10" x .75"
6120	Double decker tray base, lid, insert, 2 mats	6" x 10" x 1.5"
6510	Scope tray base, lid, and 2 silicone holders	2.6" x 16" x 1.5"
6550	Scope tray base, lid, and 2 silicone holders	6.5" x 18" x 1.5"
6560	Scope tray deep base, lid, and 2 silicone holders	6.5" x 18" x 3"
1500	Large tray w/ base, lid, & mat	10" x 15" x .75"
1520	Double decker large tray base, lid, insert, 2 mats	10" x 15" x 1.5"
1522	Large tray deep w/ base, lid, & mat	10" x 15" x 1.5"
4600	Micro tray w/ base, lid, & mat	4" x 6.5" x .75"
4700	Micro tray w/ base, lid, & mat	4" x 7.5" x .75"

*Other configurations, trays, & inserts available upon request

Stretcher Pads

All of our stretcher pads are covered in a heavy duty black Naugahyde for an extremely durable, yet comfortable feeling. Each pad comes standard with Velcro® kits and rounded corners, unless the corners are specified mitered or squared. We recommend the option of adding a layer of our [Blue Diamond®](#) gel just beneath the surface of the pad to maximize patient comfort and reduce the risk of pressure sores.



1000

Universal Stretcher Pad

Top Cut Corners w/Velcro®
24" x 76" x 3"

Available in 3" and 3.25" w/Gel
Available in 4" and 4.25" w/Gel
(See Website for Details)



1100

Standard Stretcher Pad

Rounded, Mitered or Square Corners
24"-26" x 74"-77"

Available in 3" and 3.25" w/Gel
Available in 4" and 4.25" w/Gel
(See Website for Details)




1200

Extra Wide Stretcher Pad

Rounded, Mitered or Square Corners
27"- 31" x 77"+

Available in 3" and 3.25" w/Gel
Available in 4" and 4.25" w/Gel
(See Website for Details)

Seizure/Side Rail Bumper Pads



Seizure/Side Rail Bumper Pads

Available in 3 Sizes:

Item#3800-3807: 14" x 36" x 1.5"

Item#3800-3809: 14" x 48" x 1.5"

Item#3800-3811: 14" x 72" x 1.5"

Glove Box Holders



2327M2

Glove 2-Box Holder

Powder Coated, Wall Mounted
Accessible Left, Right or Top
(Gloves Sold Separately)



2327M3

Glove 3-Box Holder

Powder Coated, Wall Mounted
Accessible Left, Right or Top
(Gloves Sold Separately)



2327M4

Glove 4-Box Holder

Powder Coated, Wall Mounted
Accessible Left, Right or Top
(Gloves Sold Separately)

Warranty Information

We back our **Blue Diamond®** Gel product line which includes Gel Bean Bag positioners with a 2 year manufacturer's defect warranty. All other O.R. positioning products shall be deemed free from any defects in material and workmanship for 1 year. What is not covered in our expressed warranty are product punctures, tears caused by misuse, neglect, and abuse.

For any additional warranty information please contact our customer service department at 800-804-0333 or your local sales representative.

Product Return Policies & Procedures

All returns will be handled through our Customer Service Department by calling 800-804-0333. If your product was received defective or was not as you ordered, please contact us immediately.

Returns will not be accepted without a pre-authorized Returns Goods Authorization Number (RGA#). Items must be returned in their original packaging and cannot have been used in an actual procedure.

To initiate a return for credit, the following criteria must be met:

- 1) Authorized RGA# must be provided to you by David Scott Company.
- 2) Product must be returned in new and reusable condition.
- 3) Product must be returned in its original packaging
- 4) Product must be returned to David Scott Company within 15 days of receipt of product.

***Custom products cannot be returned. This includes O.R. Foam Table Pads, Foam Armboard Pads, Light Cloud™ Positioners, Disposable Products, Cart Covers, and Anti-fatigue Comfort Mats.

All returns are subject to a 15% restocking fee.

Note: Refunds are for the value of the product being returned. Shipping will not be credited if the wrong product was ordered and the return shipping of said item is at your expense. Please allow up to 3 weeks for the credit to be processed.

***While we do our best to represent in this catalog product description and dimensions, David Scott Company cannot be held liable for mistakes in this catalog. Products may vary from image depictions. Please contact us to confirm all product information.

Copyright Information

Copyright 2020, David Scott Company. **Blue Diamond®** and **LightCoud™** are registered trademarks of David Scott Company. **Allen®** is a registered trademark of Allen Medical Systems, **Velcro®** is a registered trademark of Velcro Industries BV., **OSI®** and **Schlein®** are registered trademarks of Mizuho OSI, and **DeMayo®** is a registered trademark of IMP. All other trademarks are the property of their respective companies.

Product Listings

A

Accessory Cart	23
Adjustable Foot Board	20
Adjustable Stools	34
Add-A-Side-Rail	20
Anesthesia Screens	26
Flexible, Adjustable, Malleable, Table Mount	26
Anesthesia Temperature Dots	26
Anti-Fatigue Comfort Mats	33
Anti-Fog Solution	26
Armboards	
Bariatric Armboard	15
Carbon Fiber Armboard	15
C-Thru Transparent Armboard	15
Gearless Armboard	15
Kross Armboard	15
Lollipop Armboard	15
Lateral Armboard	15
Multi-Task Armboard	15
Ratcheting Armboard	15
Armboard Pads	
Foam	14
Gel	6, 14
Reusable	14
SheerCare, UltraCare & Combo	14
Arm & Hand Tables	16
Basic Carpa Table	16
Basic Hourglass Table	16
Basic Rectangle Table	16
Carbon Fiber Table	16
C-Thru Transparent Table	16
Mobile Base Table	16
Under Mount Table	16
Arm Positioners	15, 16
Arm Shields	15
Angled	15
Sled Style	15
Arm Supports	6, 14, 15, 16, 24
Arthroscopic Stress Post	27
Arthroscopic Well Leg Holder	27

B

Bariatric	28
Armboard	15, 28
Bean Bag Positioners	11, 28
Lift-Assisted Stirrups	22
Mayo Stand	28, 30
Nissen Thigh Straps	19, 28
Restraint	19
Table Width Extender	20
Beach Chairs	24
Bean Bag Positioners	
Gel & Traditional	10, 11
Accessories	11

Blue Diamond® Gel Pads	4-9
Armboard	6
Arm & Elbow	6
Arm Shield/Toboggan	6
Bean Bag	10, 11
Chest Rolls, Adult & Pediatric	7
Contoured Positioners	6
Foot, Heel, Knee & Stirrup	5
Head Positioners	4
O.R. Table Pads	8
Repair Patch	25
Specialty Pads	9
Boot Pads	
Clamshell Boot Pads	23
Foot Hugger Boot Pads	23
Gel Liner Boot Pads	5, 23
Braces	24
Arthroscopic Stress Post	27
Lateral Braces	24
Restraint Straps	19
Brachytherapy Drapes & Frame for Urocatch	21

C

Candy Cane	22
Gel Ankle Protectors	5, 22, 23
Gel Sheaths	5, 22, 23
Hooks & Straps	23
Stirrups	22
Carts	23
Accessory	23, 25
Stirrup	23
Cart Covers	35
C-Flex® & Accessories	29
Clamp, Side Rail Sockets	20

D

DaVinci Gel Bump	19
Deluxe Arthroscopic Leg Holder	27
Disposables	29
Hand, Head, Heel	14, 29
Knee, Leg	29
Spine	29
Steep Trendelenberg	29
Temperature Dots, Anesthesia	26, 29
Velcro®	25

E

Extender, Side Rail	20
Extremity Positioner	27

F

Flexible Anesthesia Screen	26
Foam Positioners	

Armboard Pads	14
Extremity Positioning Pad	27
Frogger Pad	27
Light Cloud™ Resposable	14
O.R. Table Pads	12, 13
Well-Leg Holder	27
Foot Board	20
Foot Stools	33
Fromm Knee Triangles	27

G

Gel Pads & Positioners	4-9
Armboard	6
Arm & Elbow	6
Arm Shield/Toboggan	6
Bean Bag	10, 11
Chest Rolls, Adult & Pediatric	7
Contoured Positioners	6
Foot, Heel, Knee & Stirrup	5
Head Positioners	4
O.R. Table Pads	8
Repair Patch	25
Specialty Pads	9

H

Hampers	32
Folding, Round, Square	32
Hand Positioner	16, 29
Hand Pump, Bean Bag	11
Head Block Foam	29
Head Positioners	
Foam	14, 24, 29
Gel	4
Hermanator Sterile Pads	27, 29
Hip Positioners	24
Stulberg Positioner	24
Wixson Hip Positioner	24

I

Infant Stirrups	22
Instrument Tables	31
Irrigation Tower	32
IV Poles	32
Chrome	32
Easy Lift Hand Crank	32
Irrigation	32
Power Strip	32
Stainless Steel	32

J

Jackson Table Covers	29
----------------------	----

Product Listings

K

Kick Bucket & Frame	32
Knee Positioners	27
Knee Crutch & Accessories	27
Knee Replacement System (TKR)	27
Knee Stabilizer, Total (TKS)	27
Kross Armboard & Cloth	15

L

Large Patient, See Bariatric	28
Lateral Brace System	24
Lateral Positioners	24
Bean Bags	10, 11
Lateral Brace System	24
Peg Board	17
Stulberg Positioners	24, 27
Wixson Hip Positioner	24
Leg Holders	27
Arthroscopic	27
DaVinci Bump	19
Leg Prepper, Single Loop	27
Nissen Straps	19
Restraint Straps	19
Leg Prepper	27
Light Cloud™ Resposable Positioners	14
Head Donuts, Heel, Prone	14

M

Mat, Comfort	33
Mayo Stands	30
Multi-Task Armboard	15

N

Nerve Stimulators	26
Alligator Leads	26
Diagnostic Ball Probe	26
Protective Covers	26
Neonatal Gel Positioners	4, 7
Nested Tables	31
Nissen Thigh Straps	19
Nissen Thigh Straps, Bariatric	19

O

Offset Accessory Rail	20
OR Table Pads/Surfaces	12, 13
SheerCare, UltraCare, Combo	12, 13

P

Pads	
Armboard Pads	14
Clam Shell Boot Pads	23
Frogger Pad	27

Gel Pad	4-9
OR Table Pads	12, 13
Standard OR Table Pads	12, 13
Stirrup Pads	5, 23
Table Pads	8, 12, 13
Patient Shifter Boards	18
Patient Transfer Boards	18
Pediatric	
Gel Positioners	4, 7
Stirrups	22
Peg Boards & Accessories	17
Picket Fence Leg Holders	27
Plastic Sterilization Trays	35
Prepper Device	27
Prone Positioners	
Arm Support	6, 15
Gel Table Pads	8
Head Positioners	4, 14, 29

R

Rail Clamp	20
Repair Patches, Gel & Bean Bag	25
Ratcheting Armboards	15
Resposables, Light Cloud™	14
Restraining Gel Bump	19
Restraint Straps & Accessories	19
Robb Leg/Knee Positioner	27
Roller Boards & Accessories	18

S

Samarit Roller Boards	18
Security Straps	19
Seizure Pads	36
Shifter Boards	18
Short Accessory Rail	20
Shoulder Positioners	24
Lateral Braces	24
Shoulder Brace Pads	9
Shoulder Beach Chair	24
Shoulder Cutout Bean Bag	10, 11
Shoulder Chair & Accessories	24
Shoulder Holder	24
Shoulder Supports	24
Side Rail Bumper Pads	36
Side Rail Extender	20
Side Rail Sockets/Clamps	20
Simple Clamp	20
Single Loop Leg Prepper	27
Sled-Style Arm Shields	15
Spine	25
Disposables	29
Schure Spine Frame	25
Wilson Frame	25

Stands	
Hamper	32
Mayo	30
Solution	30
Step Stools	33
Sterilization Trays, Plastic	35
Steep Trendelenberg Positioner	29
Stirrups	
Bariatric	22
Basic	22
Candy Cane	22
Cart Holder	23
Gel Boot Liner Pads	5, 23
Lift-Assisted	21, 22
Pediatric	22
Replacement Pads, Foam	23
Stirrup Carts	23
Stools, Adjustable	34
Straps	19
Stretcher Pads	36
Stress Post, Arthroscopic	27
Stulberg Positioners	24, 27
Suction, Bean Bag	10, 11

T

Table	
Instrument, Nested, Utility	31
Table Width Extenders	20
Temperature Dots, Anesthesia	26
Thigh Straps, Nissen	19
Toboggan Arm Shields	15
Total Knee Stabilizer	27
Transfer Boards	18

U

Universal Accessory Clamp	20
Urology Urocatch Drapes & Frame	21
Utility Tables	31

V

Vacuum Pac Positioner	10, 11
Valve, Bean Bag Replaceable	11
Vacuum Pac Hand Pump	11
Velcro®	19, 25

W

Well-Leg Holder, Arthroscopic	27
Well-Leg Holder, Foam	27
Wilson Frame, Radiolucent	25
Wixson Hip Positioner	24

DAVID SCOTT
COMPANY

2 Master Drive • Franklin MA 02038
508.875.3333 • 800.804.0333 • Fax: 508.875.3375

Blue  Diamond®

alpinemédica ®

alpinemédica es: Jesús Del Pino Favila PIFJ731124UW5
Oficinas: Jaime Gomez No. 2802 Int. 2 Col. Paseos Del Sol
Zapopan, Jalisco, México, C.P. 45079
Tels: (+52) 33-15-78-28-02 / 33-15-68-16-02
www.alpinemédica.com.mx