



## Bed Techs Inc.

Specializes in Reconditioned Hill Rom® Beds & Stretchers

### Hill Rom® Advanta

Hill Rom® Beds for Safety, Quality, Comfort. Great Products—Intelligent, Well-Designed, Clinically Effective and Easy to Use.

#### Our Process Includes

- Professionally refurbished medical beds to manufacturer's specifications
- JCAHO required documentation
- New paint & labels
- Fully function tested
- Reconditioned PC boards and side rail switches
- Reconditioned side rails
- Reconditioned head/foot boards
- Fresh hydraulic fluid
- Rebuilt casters with new components and wheels
- Delivery/Installation performed by our technicians (Domestic only)
- 1 year warranty (Domestic only)
- Quality & Price

#### Specifications

- Side communications (Nurse call, lighting, entertainment, volume and channel controls)
- Bed exit
- Zone air mattress
- Tente casters
- Scale
- Bumpers
- Automatic night light
- Instant CPR release
- Height in high position 30"
- Height in low position 17"
- Sleep surface width 36"
- Sleep surface length 82"
- Overall length (w/out bumpers) 36" x 91"

#### Options

- Service Contracts for preventative maintenance
- 220-volt transformers

